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Estimation of the total hardness of a water sample using edta

Theory


Ethylenediaminetetraacetic acid (edta) is a reagent that forms edta-metal complexes with many metal ions (**but not with alkali metal ions such as Na⁺ and K⁺**). In alkaline conditions (pH>9) it forms **stable complexes** with the alkaline earth metal ions **Ca²⁺ and Mg²⁺**. The edta reagent can be used to measure the total quantity of dissolved

Ca^{2+} and Mg^{2+} ions in a water sample. Thus the total hardness of a water sample can be estimated by titration with a standard solution of edta.

Suitable conditions for the titration are achieved by the addition of a buffer solution of pH 10. The buffer solution stabilises the pH at 10. There are H^+ ions produced as the reaction proceeds, and without the buffer solution the pH would decrease.

The edta reagent cannot under these conditions distinguish between the hardness caused by Ca^{2+} and Mg^{2+} , or (directly) between temporary and permanent hardness. Therefore the results of this experiment are usually expressed in terms of the quantity of insoluble CaCO_3 that would have to be converted into soluble salts to give the same total number of moles of dissolved Ca^{2+} and Mg^{2+} ions. This enables the total hardness of water from different sources to be compared easily.

Because it is a primary standard, and is also more soluble in water, the disodium salt of edta is more commonly used as the reagent rather than edta itself. If $\text{Na}_2\text{H}_2\text{Y}$ represents this salt, it ionises in aqueous solution to H_2Y^{2-} , which complexes in a 1:1 ratio with either Ca^{2+} or Mg^{2+} ions (which are represented as M^{2+}). The reaction can be represented as follows:

The indicator Eriochrome Black T is used to detect the end point. This is an indicator that has a different colour when complexed to metal ions than when it is a free indicator. The reaction between the red indicator-metal complex and the edta reagent at the end point can be represented as follows:

Wine red blue

Chemicals and Apparatus

Hard water sample

0.01 M edta solution

Buffer solution, pH 10

Eriochrome Black T indicator ("diluted" with sodium chloride)

Deionised (or distilled) water


Pipette (25 cm³)

Pipette filler
 Burette (50 cm³)
 Conical flask (250 cm³)
 Filter funnel
 Beakers (250 cm³)
 Filter paper or white tile
 White card
 Wash bottle
 Graduated cylinder (10 cm³)
 Spatula for solid indicator
 Retort stand
 Boss head
 Clamp
 Safety glasses

Procedure

Wear your safety glasses.

1. Wash the pipette, burette and conical flask with deionised water. Rinse the burette with the edta solution and the pipette with the hard water.
2. Using the funnel, fill the burette with the edta solution. Open the tap briefly to fill the part below the tap. Remove the funnel. Adjust the level of the solution to the zero mark. Make sure that the burette is vertical.
3. Use the pipette to transfer 50 cm³ of the hard water sample to the conical flask.
 Add 2-3 cm³ of the buffer (pH 10) solution (measured out using the graduated cylinder).
4. Add 0.03 g of the solid indicator to the contents of the flask in the following manner:
 Add gradually to the flask, swirling after each addition. A deep wine red colour is obtained.
5. Carry out one 'rough' titration to find the approximate end point, followed by a number of accurate titrations until two titres agree to

within 0.1 cm³. At the end point, the colour should be dark blue, with no tinge of wine-red colour.

6. From the data, calculate the total hardness of the water sample.

Table of Results

Volume of hard water sample	=	cm ³
Molarity of edta solution	=	M
Rough titre	=	cm ³
Second titre	=	cm ³
Third titre	=	cm ³
Average of accurate titres	=	cm ³
Total hardness	=	mol/l Ca ²⁺
Total hardness	=	g/l CaCO ₃
Total hardness	=	ppm CaCO ₃

Questions relating to the experiment

1. Why is it important that the reaction between the edta and the metal ions in solution (i) is rapid and (ii) goes to completion?
2. The water sample could contain metal ions other than Ca²⁺ and Mg²⁺. How would the reliability of the result be affected if this were the case? Suggest two other metal ions that could be present in the water.
3. This reagent cannot distinguish between temporary and permanent hardness. List the compounds of calcium and magnesium that cause hardness, and indicate those which cause temporary hardness.
4. Suggest a method of establishing the amount of permanent hardness in a water sample.
5. What is the function of the buffer solution?

Forrás:

https://docs.google.com/document/d/1DT3szwPMdh_s11GR0ho5zAWWwjIlBB9Civ6SP9f1u48/edit